

ITQB-Nova - Universidade Nova de Lisboa	EPP-SOP-ITQB01
	Version 1.0

EPP - Standard Operating Procedure

(only for selected experiments intended to transfer results from one lab to the other)

Title: Extraction of total RNA from *P. putida*

Total RNA extraction for downstream applications

<u>distribution list</u>			
changes to prior version:			
	name	signature	date
experimenter 1	Susana Domingues		27.04.2017

Instruction

Extraction of total RNA from *P. putida*

1 Introduction / Purpose

This protocol describes the method for total RNA extraction from *Pseudomonas putida*. RNA obtained from this method ranges from large transcripts to small RNA molecules. It is not completely DNA free, so Real-Time PCR downstream applications need further RNA sample treatment with a DNase.

2 Equipment and chemicals

2.1 Equipment

- ⇒ Centrifuge
- ⇒ Water bath
- ⇒ Tissue lyser (FastPrep)

2.2 Chemicals

Tris base
Lysozyme
Glass beads
Acid phenol (for RNA)
DNase RNase free (Turbo DNase Ambion)
Trizol
CHCl₃
Isoamyl alcohol
isopropanol
Ethanol (EtOH)

2.3 Bacterial strains

This protocol can be applied to any *P. putida* strain

3 Media and buffers

Lysis solution: Weight 20 mg of lyophilized powder lysozyme and mix in 10 ml of Tris-Cl buffer 10mM pH7.5. Aliquot the solution and store -20 °C.

4 Procedures

4.1 Cell Lysis

Resuspend the collected pellets into 350 µl of lysis solution in 2 ml tube and incubate for 5 min at 37°C.

Change the suspension to a 2ml eppendorf with screw cap containing 0.5ml glass beads (smaller ones) 0.5ml acid phenol (RNA).

Lyse the cells in the Tissue lyser (FastPrep) – speed 6, 45s (2x).

4.2 DNase Treatment

Centrifuge max speed, 10min, 4°C.

Take the aqueous phase into 2ml eppendorfs

Add 50µl DNase Buffer 10x

Add 10U DNase RNase free (Turbo DNase Ambion). Incubate at 37°C, 1H.

4.3 RNA extraction

Add 1ml Trizol + 200µl CHCl₃

Vigorous agitation 2 min room T

Centrifuge 14000 rpm, 10min, 4°C

Take the aqueous phase

Add 200µl CHCl₃:Isoamyl alcohol (24:1)

Vigorous agitation 2 min room T

Centrifuge 14000 rpm, 15min, 4°C

Take the aqueous phase

Add 500µl isopropanol

Soft vortex 5 to 10sec

10min at room T

Centrifuge 14000 rpm, 8min, 4°C

Wash with EtOH 70%

Centrifuge 14000 rpm, 5min, 4°C

Dry at 37°C ~15min

Resuspend in 100µl H₂O (let hydrate before resuspending)

5 *remarks/troubleshooting:*

6 *Biosafety*

Work in a fume hood

Wear protective gloves and lab coat.

7 *Acknowledgements*



This project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 635536.