

WCSB, Wageningen University	IBISBA-SOP-WU26
	Version 1.0

## EPP - Standard Operating Procedure

(only for selected experiments intended to transfer results from one lab to the other)

Title: *Marker-less multiple gene deletions in P. putida KT2440*

<u>distribution list</u>			
changes to prior version:			
	name	signature	date
Experimenter 1	Rita Volkers		19/3/2019

## Instruction

### *Marker-less multiple gene deletions in P. putida KT2440*

#### 1. Introduction / Purpose

In this protocol a method is described to make multiple gene-deletion in *P. putida* KT2440. The protocol is an adaptation for *P. putida* of the method described by Fehér *et al.* (for *E. coli*). The two plasmids described here, pSW(I-Sce-I) and pEMG, can be delivered into the cells by electroporation or tri-parental mating (using an *E. coli* helper strain that harbors the pRK600 plasmid).

For this deletion process the I-SceI expressing plasmid (pSW) could be introduced before (preferred option if you plan to do multiple deletions within the same strain) or after obtaining your co-integrate. The second will be described here.

**Keywords:** Gene deletion – marker-less – KT2440 – *P. putida* – knock-out

#### 2. Equipment and chemicals

##### 2.1. Equipment

- Centrifuge for eppendorf tubes
- Electroporation machine
- PCR machine

##### 2.2. Chemicals

- PCR purification kit
- PCR reagents
- 300 mM sucrose
- 1-2mM 3-methylbenzoate

##### 2.3. Bacterial strains

*Pseudomonas putida* KT2440, *E.coli* Dh5αpir

##### 2.4. Other materials

- Plasmid pSW(I-Sce-I)
- Plasmid pEMG or pGNW
- LB plates with kanamycin, x-gal, and iptg
- LB plates with kanamycin 50 µg ml<sup>-1</sup> and ampicillin 600 µg ml<sup>-1</sup>
- LB plates with ampicillin 600 µg ml<sup>-1</sup>
- LB plates with ampicillin 500 µg ml<sup>-1</sup>
- Primers outside of TS1 and TS2
- M13 primers: *fw=tgtaaaacgacggccagt* and *rev=caggaaacagctatgaccatg*

#### 3. Media and Buffers

- LB medium

#### 4. Procedures

The current example explains the procedure for the deletion of *aceEF* operon (*Pseudomonas putida* KT2440).

Cloning region of interest into the suicide plasmid pEMG or pGNW

The cloning step could be performed by 2 methods: [1] PCR regions and cloning each region separately into the plasmid or [2] by overlap PCR, which is described here.

##### Overlap PCR

To construct the corresponding recombinogenic vectors, 500 bp (TS1) upstream and 500 bp downstream (TS2) regions of *aceEF* operon should be amplified separately and then join in a single DNA segment by overlap PCR (Fig. 1)

- Do the first PCR reactions to amplify the 500 bp upstream and downstream regions using a fresh colony of KT2440 as template.
- Prepare the second overlap PCR reaction using 1µl from each of the reaction products of the first PCR as template.
- Clean up the second PCR product using a purification column, digest with the appropriate restriction enzymes, inactivate them, and ligate to the suicide vector (pEMG or pGNW).




Figure 1

##### Transformation of *E. coli* Dh5αpir

- Transform *E. coli* DH5αpir with an aliquot of the ligation mixture and plate onto LB plates with kanamycin, X-Gal and IPTG.
- Select a few white colonies, re-streak them and check for the presence of the correct plasmid.
- Use M13 primers (few bases outside of the pGNW/pEMG polylinker) to select a positive clone and send it to sequence to confirm that the targeting sequences don't have sequence errors.
- **Note:** Use *I-SceI* and the restriction profile of the constructed plasmid in the case that you use M13 primers for selection and sequencing. M13 bind just a few bases far from *I-SceI* sites of pEMG, so the sequencing cannot give you reliable results for that region

Cointegration of pEMG-aceEF(TS1-2)

*Preparation of P. putida KT2440 electrocompetent cells and introduction of pEMG-aceEF(TS1-2) or pGNW-aceEF(TS1-TS2)*

- Grow overnight a 20 ml culture of KT2440 in LB medium at 30 °C with shaking.
- Centrifuge at room temperature for 10 minutes.
- Remove supernatant and resuspend the pellet gently in 10 ml of 300 mM sucrose and centrifuge as in step 2.
- Discard supernatant, add 1 ml of 300 mM sucrose, resuspend the pellet and transfer it to an eppendorf tube.
- Centrifuge at 12000 rpm for 2 minutes at room temperature.
- Dispose of supernatant and add 500 µl of 300 mM sucrose, resuspend and distribute in 100 µl aliquots.
- Add 300-400 ng of your construct plasmid to a 100 µl aliquot of KT2440, mix gently and transfer to a 2 mm gap width electroporation cuvette and proceed to electroporate.
- Add 1 ml of LB or SOC and incubate 2 hours at 30 °C with shaking.
- Plate everything onto LB+kanamycin 50 µg ml<sup>-1</sup> and incubate the plates for 2 days. (After 2 days there are a few big colonies which probably contain 2 population, the Km<sup>r</sup> and Km<sup>s</sup>)
- Streak a few big colonies and use colony PCR to check for co-integration events.

**Note:** For the selection of the right co-integrates, design and use primers which bind outside of the TS1-gene-TS2 region. In this way you are looking just for a single large PCR product instead of selecting with TS1-Fwr, TS2-Rvr. In the case of aceEF the desired product had the size of 10kb (Fig. 2 and Fig. 3).




Figure 2




Figure 3

### Delivery of pSW (I-SceI)

- Prepare electrocompetent cells of your co-integrate as described before.
- The pSW(I-SceI) plasmid described here bears the ampicillin resistance gene.
- Add 50 ng of pSW(I-SceI) plasmid to a 100  $\mu$ l aliquot of your co-integrate, mix gently and transfer to a 2 mm gap width electroporation cuvette and proceed to electroporate.
- Add 1 ml of LB or SOC and incubate 1 hour at 30 °C with shaking.
- Plate several dilutions onto LB+Km50+ampicillin 600  $\mu$ g ml<sup>-1</sup> (KT2440 is naturally resistant to ampicillin but such concentration allows plasmid selection).
- Streak a few colonies and check the presence of the plasmid by miniprep or PCR (not essential step).

### Induction of the I-SceI enzyme

- Start by diluting overnight cells to an OD600 of 0.1 and induce with 1-2mM 3-methylbenzoate for 6-7 h.
- Plate a few dilutions onto LB or LB plus Amp600 (if you wish to maintain the pSW-I plasmid for further deletions).
- Restreak colonies in LB and LB+Km50 to check the loss of the co-integrated plasmid.
- Select a few kanamycin sensitive clones and differentiate by colony PCR between WT and deleted clones (the proportion is usually close to 50%). For the selection of the knockout colonies use the primers outside of TS1 and TS2. The knockout colonies should give a size product of the size of TS1-TS2 (col 5-6) and the wild type the normal control size (col 1, 2, 3, 4) (Fig. 4).



Figure 4

#### Curation of pSW(I-SceI)

- Pick a clone with deletion and do several passages in LB without antibiotics. In order to do that grow overnight in LB at 30 °C with shaking, then transfer to a new tube with fresh LB and repeat the process for 3 days.
- Plate onto LB appropriate dilutions to obtain separate colonies. Restreak several colonies on LB and LB+amp500.
- Select amp sensitive colonies and double check the plasmid curation by colony PCR using the primer pair described before.

#### **5. Remarks / troubleshooting**

- Use I-SceI and the restriction profile of the constructed plasmid in the case that you use M13 primers for selection and sequencing. M13 bind just a few bases far from I-SceI sites of pEMG, so the sequencing cannot give you reliable results for that region
- For the selection of the right co-integrates, design and use primers which bind outside of the TS1-gene-TS2 region. In this way you are looking just for a single large PCR product instead of selecting with TS1-Fwr, TS2-Rvr. In the case of aceEF the desired product had the size of 10kb.

#### **6. References**

Fehér T., Karcagi I., Gyórfy Z., Umenhoffer K., Csörgő B., Pósfai G. (2008) *Scarless Engineering of the Escherichia coli Genome*. In: Osterman A.L., Gerdes S.Y. (eds) *Microbial Gene Essentiality: Protocols and Bioinformatics*. Methods in Molecular Biology™, vol 416. Humana Press

#### **7. Acknowledgements**



IBISBA (grant agreement No 730976) and EmPowerPutida (grant agreement No 635536) projects have received funding from the European Union's Horizon 2020 research and innovation programme.